

CITEST COVID-19 Antigen Rapid Test (swab/nasal swab) (1 pc)



Product name	CITEST COVID-19 Antigen Rapid Test (swab/nasal swab) (1 pc)
Article number	04-024
Product description	<p>CE-certified as self-test for self-testing (CE1434) BfArm listing under AT1350/21 Overall sensitivity: 100% / Specificity: 99.9% Test is performed via nasal swab Easy to use, result in 15 minutes All test components are included Test winner among the self-tests according to Stiftung Warentest (German consumer organization) Best rapid antigen test according to PEI evaluation (Paul Ehrlich Institute)</p>
Manufactures	CITEST® Diagnostics inc.
Country of origin	Canada
Components	<p>1x Test cassette 1x Extraction buffer 1x Sterile swab 1x Biosafety bag 1x Tube holder 1x Package insert</p>
Certificates	<p>CE1413 EC Certificate No. 1434-IVDD-054/2022</p>
Packaging information	<p>480 pieces/carton 5.760 pieces/pallet</p>
Shipping information	Cartons per pallet: 12