

COVID-19 Antigen Rapid Test (Nasal swab) Joinstar® (1 ps)



Product name	COVID-19 Antigen Rapid Test (Nasal swab) Joinstar® (1 ps)
Article number	04-021
Product description	CE-certified as self-test for self-application (CE1011) Anterior-nasal (short nose!) Sensitivity: 100% / Specificity: 98.1% All test components are included Easy to use, result in 15 minutes Evaluated by Paul-Ehrlich-Institute (PEI)
Manufactures	Joinstar Biomedical Technology Co., Ltd.
Country of origin	China
Components	1x Test cassette 1x Extraction buffer 1x Sterile swab 1x Biosafety bag 1x Tube holder 1x Package insert
Certificates	CE1011
Packaging information	400 pieces/carton 6.400 pieces/pallet
Shipping information	Cartons per pallet: 16