

FFFP3 respirator with breathing valve Opharm® VE10 (black)



Product name	FFP3 respirator with breathing valve Opharm® VE10 (black)
Article number	01-0201.1
Product description	5-layer FFP3 filtering half mask with breathing valve in black Unisex size Flexible metal strip at nose level BFE filtration > 99 % Personal protective equipment according to EN149:2001 + A1:2009 Packed 10 masks/pack Made in EU
Manufactures	Opharm® S.A.
Country of origin	Poland
Components	10 FFP3 respirator with breathing valve
Certificates	ISO 9011:2015 EN 149: 2001 + A1: 2009 CE certified
Packaging information	240 pieces/carton 8.640 pieces/pallet
Shipping information	Carton size: 30cm x 40cm x 40cm Cartons per pallet: 36