

FFP2 NR mask Opharm® VE5 (black)



Product name	FFP2 NR mask Opharm® VE5 (black)
Article number	01-0010.1-VE5.4
Product description	5-layer FFP2 filter half mask in black color Black color Unisex size Flexible metal strip at nose level BFE filtration > 95 % Personal protective equipment according to EN149:2001 + A1:2009 Packed 5 masks/pack Made in EU
Manufactures	Opharm® S.A.
Country of origin	Poland
Components	5 FFP2 NR masks in black
Certificates	ISO 9011:2015 EN 149: 2001 + A1: 2009 CE certified
Packaging information	450 pieces/carton 16.200 pieces/pallet
Shipping information	Carton size: 30cm x 40cm x 40cm Cartons per pallet: 36