

Nitrile disposable gloves Opharm® (100s) - Black



Product name	Nitrile disposable gloves Opharm® (100s) - Black
Article number	03-073-0
Product description	Medical, diagnostic disposable protective gloves, non-sterile. Textured fingertips provide secure grip in a variety of situations, increase friction and grip precision. Latex-/Powder-free 4g nitrile gloves for work in a variety of settings. Provides protection against viruses, chemicals and microorganisms.
Manufactures	Opharm S.A.
Country of origin	Poland
Components	100 disposable nitrile gloves
Certificates	EN ISO 374-5:2016, EN ISO 374-1:2016 / Type B, EN 420:2003+A1:2009, EN 455 1, 2, 3, 4
Packaging information	100 pcs/box 10 boxes/carton