

## TypIIIR medical mask Opharm® VE50 (black)



<b>Product name</b>	TypIIIR medical mask Opharm® VE50 (black)
<b>Article number</b>	01-0201.1
<b>Product description</b>	3-layer medical mask typIIIR Unisex size Flexible metal strip at nose level BFE filtration > 98 % Certification CE & EN 14683:2019 Packed 50 masks/pack Made in EU
<b>Manufactures</b>	Opharm® S.A.
<b>Country of origin</b>	Poland
<b>Components</b>	50 typIIIR medical mask
<b>Certificates</b>	CE certified EN 14683:2019 EN 1041:2008+A1:2013 EN ISO 14971:2012 ISO 9001 ISO 13485
<b>Packaging information</b>	1.200 pieces/carton 43.200 pieces/pallet
<b>Shipping information</b>	Carton size: 30cm x 40cm x 40cm Cartons per pallet: 36