

## Solmira® Combotest 4 in 1



<b>Product name</b>	Solmira® Combotest 4 in 1
<b>Article number</b>	04-027
<b>Product description</b>	Detects RS virus, Influenza A/B, Covid-19 Self-test Combo 4 in 1 Result in 15 minutes Easy to use 18 months shelf life Sensitivity: 92.93% Specificity: 100%
<b>Manufactures</b>	Shenzhen Microprofit Biotech Co., Ltd.
<b>Country of origin</b>	China
<b>Components</b>	1 sterile swab 1 extraction tube with buffer solution integrated incl. drop caps 1 test cassette 3 compartments 1 EN instruction manual
<b>Certificates</b>	CE Certification 2934
<b>Packaging information</b>	500 pcs./carton