

Latex disposable gloves Top Glove® (100s) - Natural



Product name	Latex disposable gloves Top Glove® (100s) - Natural
Article number	03-078
Product description	Latex gloves with medical suitability, Good abrasion and tear resistance, Ambidextrous examination gloves, Suitable for medical use, Approved for food contact, Powder-free, Non-sterile, AQL 1.5
Manufactures	Top Glove
Country of origin	Malaysia
Components	100 disposable latex gloves
Certificates	EN ISO 374-5:2016, EN ISO 374-1:2016 / Type B, EN 420:2003+A1:2009, EN 455 1, 2, 3, 4
Packaging information	100 pcs/box 10 boxes/carton